Clinical specimens for testing at the NCSLPH should be shipped within 24 hours of collection.

When shipping must be delayed (e.g., inclement weather, courier closures), guidelines are given below on proper storage and shipping of common specimen types. Please see our <u>SCOPE</u> for detailed guidance.

	If shipping immediately	If shipping is delayed	
COVID-19 + Flu PCR	≤ 72 hours	🎇 DRY ICE REQUIRED 🎇	
Legionella Clinical**	≤ 3 days	🎇 DRY ICE REQUIRED 🎇	
CHIKV, EEE, LACV, SLE, WEE, WNV	≤ 2 days	🎇 DRY ICE REQUIRED 🇱	
Dengue (DENV)	≤ 5 days	🎇 DRY ICE REQUIRED 🇱	
Hepatitis A (HAV)	≤ 5 days	🎇 DRY ICE REQUIRED 🇱	
Hepatitis B (HBV)	≤ 4 days	🎇 DRY ICE REQUIRED 🎇	
Hepatitis C (HCV)	≤ 5 days	🎇 DRY ICE REQUIRED 🎇	
HIV	≤ 5 days	🎇 DRY ICE REQUIRED 🎇	
HIV/HCV	≤ 5 days	🎇 DRY ICE REQUIRED 🎇	
Measles/Mumps PCR*	≤3 days	🎇 DRY ICE REQUIRED 🎇	
Measles Serology*	≤3 days	🎇 DRY ICE REQUIRED 🎇	3
Rickettsia + Ehrlichia	≤3 days	🎇 DRY ICE REQUIRED 🎇	7145
Rubella IgG	≤ 5 days	🎇 DRY ICE REQUIRED 🎇	
Syphilis RPR	≤5 days	🎇 DRY ICE REQUIRED 🎇	
Syphilis TP	≤ 5 days	🎎 DRY ICE REQUIRED 🎎	
Bordetella Culture**	≤3 days		
Bordetella PCR**	≤7 days		**
HSV/VZV NAAT	≤ 5 days		
Hemoglobinopathy Whole Blood	≤6 days		
Legionella Culture**	Ship as soon as possible		
Blood Lead	≤ 43 days, including prenatal		
Micro Reference Cultures	Ship ambient as soon as possible		
NBS Dried Blood Spots (DBS)	\leq 14 days, avoid excessive heat and/or humidity		
Hemoglobinopathy Screening DBS	≤ 14 days		
Neisseria Culture	Must be shipped ambient as soon as possible, see SCOPE for N. gonorrhoeae subculture details		
BTEP* and CTAT*	Call 919-807-8600 (BTEP) or 919-602-2481 (CTAT) prior to any submission for guidance on collection, labeling, packaging, and shipment		
Chlamydia/Gonorrhea TMA	May be received cold or ambient, ≤ 28 days		
Enteric Bacteriology	Varies by specimen type, see SCOPE for details		
Mycobacteriology (TB)	Varies by specimen type, see SCOPE for details		
Mycology	Varies by specimen type, see SCOPE for details		
Parasitology	Varies by specimen type, see SCOPE for details		
Rabies	Ship cold ASAP, freezing not recommended		
Viral Culture	Varies by specimen type, se	e SCOPE for details	

LHDs follow NCSLPH DOA Medical Courier Packaging & Shipping Guidance.

Refrigerate specimens 2-8°C post collection. Specimens must be received cold (2-8°C) on frozen ice packs within the number of hours/days as indicated in the chart from time of collection.

- Freeze specimens at ≤ -70°C post collection. Ship overnight via commercial courier on dry ice.
 - Freeze specimens at ≤ -20°C post collection. Ship overnight via commercial courier on dry ice.
 - ** Specimens must be received cold (1-8°C) on frozen ice packs * Submission requires prior approval

Store and ship at ambient temperature. Specimens must be received within the number of hours/days as indicated in the chart from time of collection.

Unique requirements, refer to SCOPE.



NC Department of Health and Human Services Division of Public Health North Carolina State Laboratory of Public Health

https://slph.dph.ncdhhs.gov/

NCDHHS is an equal opportunity employer and provider 7/18/2023

NCSLPH DOA Medical Courier Packing and Shipping Guidance