# COVID-19 & Influenza Specimen Collection and Shipment

## North Carolina State Laboratory of Public Health

All COVID-19 & Influenza specimens submitted to the NCSLPH must meet the testing criteria on our website and be <u>collected safely</u>. The following guidance applies only to the NCSLPH.

◆ Failure to follow submission guidelines may result in specimen rejection.

#### **Specimen Collection**

The following are acceptable specimen types:

Nasopharyngeal (NP) swab Oropharyngeal (OP) swab Nasal swab Mid-turbinate swab

Use a synthetic tipped, sterile swab of appropriate size with a plastic or metal shaft. Do not use a calcium alginate swab or one with a wood shaft.

The following are acceptable media types:

VTM (i.e., Remel M4RT) MTM (i.e., PrimeStore)

Label the vial containing the swab completely:

Swab Type
Patient Name
Date of Birth
Or any 2 unique identifiers
Date of Collection

Store properly until shipment:

Refrigerate at 2-8° C for up to 72 hours Freeze ≤ -70° C for longer storage

 Completely fill out the COVID-19 submittal form online in our ETOR portal.

#### Specimen Shipment

 All specimens shipped commercially should be packaged <u>UN3373</u> Category B, Biological Substance



If using a private courier or LHD medical courier, follow pack and ship requirements of the courier.

- ★ Keep cold. Ship refrigerated specimens on frozen cold packs and any ≤ -70° C frozen specimens on dry ice (the latter following IATA regulations). If shipping using a medical courier, no cold packs are required but specimens must be refrigerated.
- Label the package completely:

Attention: Virology/Serology Unit
North Carolina State Laboratory of Public Health
4312 District Drive
Raleigh, NC 27607-5490

 Specimen deliveries are only received at the State Lab loading dock from:

8 a.m. - 5 p.m., Monday - Friday 8 a.m. - noon, Saturday

### **Result Reporting**

Results are posted electronically to our <u>CELR</u> and <u>ETOR</u> systems, tied to the submitter's EIN number. Visit our <u>website</u> for account setup and tutorials.

