



NORTH CAROLINA STATE LAB OF PUBLIC HEALTH  
 CHEMICAL TERRORISM SERUM SPECIMEN COLLECTION  
 AND SHIPPING MANIFEST

PLACE A  $\checkmark$  IN EACH BOX FOR SAMPLES SHIPPED-PLACE AN X IN EACH BOX FOR SAMPLES NOT SHIPPED

PLEASE INDICATE THE SIZE OF THE TUBE COLLECTED (7 OR 10 mL), AND THE DATE/TIME THAT THE SPECIMEN WAS COLLECTED IN THE COMMENTS

SERUM SPECIMENS SHOULD BE STORED FROZEN AS SOON AS POSSIBLE POST COLLECTION AND THROUGHOUT THE SHIPMENT PROCESS

RT = RED-TOP OT= ORANGE-TOP

Patient/Victim ID Label	RT 1	RT 2	OT 1	OT 2	Comments:
					_____ _____ _____
					_____ _____ _____
					_____ _____ _____
					_____ _____ _____
					_____ _____ _____

Note: Please include 2 empty red-top tubes and 2 empty orange-top tubes from each lot number collected for background measurement.

Vacutainer® type may be altered based on availability in an emergency response situation.

Packed by(print):	Signature:
-------------------	------------